

APE848Hu01 100µg

Active Diacylglycerol Kinase Epsilon (DGKe)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1th Edition (Apr. 2016)

#### [PROPERTIES]

Source: Prokaryotic expression.

Host: E. coli

Residues: Leu243~Gly502 Tags: N-terminal His-tag

**Purity: >95%** 

Buffer Formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.05% sarcosyl

and 5% trehalose.

**Applications:** Cell culture; Activity Assays; In vivo assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 7.7

Predicted Molecular Mass: 32.5kDa

Accurate Molecular Mass: 33kDa as determined by SDS-PAGE reducing conditions.

#### [USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.



**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

### [SEQUENCE]

LLNPVQVF

DVTKTPPIKA LQLCTLLPYY SARVLVCGGD GTVGWVLDAV DDMKIKGQEK YIPQVAVLPL GTGNDLSNTL GWGTGYAGEI PVAQVLRNVM EADGIKLDRW KVQVTNKGYY NLRKPKEFTM NNYFSVGPDA LMALNFHAHR EKAPSLFSSR ILNKAVYLFY GTKDCLVQEC KDLNKKVELE LDGERVALPS LEGIIVLNIG YWGGGCRLWE GMGDETYPLA RHDDGLLEVV GVYGSFHCAQ IQVKLANPFR IG

## [ACTIVITY]

DGKe (Diacylglycerol kinase epsilon) is an enzyme that mainly involved in the regeneration of phosphatidylinositol (PI) from diacylglycerol in the PI-cycle during cell signal transduction. PDHA1 (Pyruvate dehydrogenase E1 component subunit alpha) has been proven to be an interactor of DGKe through Affinity Capture-MS. Thus a binding ELISA assay was conducted to detect the interaction of recombinant human DGKe and recombinant human PDHA. Briefly, DGKe were diluted serially in PBS, with 0.01%BSA (pH 7.4). Duplicate samples of 100uL DGKe were then transferred to PDHA-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-DGKe pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50µL stop solution to the wells and read at 450nm immediately. The binding activity of DGKe and PDHA was shown in Figure 1, and this effect was in a dose dependent manner.




Figure 1. The binding activity of DGKe with PDHA.

# [ IDENTIFICATION ]




Figure 2. SDS-PAGE

Sample: Active recombinant DGKe, Human




Figure 3. Western Blot

Sample: Recombinant DGKe, Human;

Antibody: Rabbit Anti-Human DGKe Ab (PAE848Hu01)